

## MICRONIC | SAMPLE HANDLING | CAPPING & DECAPPING

### MANUAL CAPMAT SEALER

DO YOU REQUIRE A SIMPLE AND COMPLETE SOLUTION TO CLOSE YOUR SAMPLES WITH TPE PLUG STYLE CAPS? THEN, THE MANUAL CAPMAT SEALER IS THE IDEAL MANUAL CAPPING SOLUTION FOR YOU. SAMPLE INTEGRITY IS A KEY REQUIREMENT FOR THE LONG-TERM STORAGE OF VALUABLE SAMPLES.

The Micronic Manual Capmat Sealer guarantees an equal sealing of plug style tubes in a 96-well configuration. This manual device requires no installation and is simple to operate.

The Manual Capmat Sealer can be used to seal any TPE cap on the Micronic Roborack-96, Lobarack-96 and Comorack-96. Using racks with an SBS-footprint minimises the storage space of your samples.

Standard for 1.40ml tubes in all configurations like: non-coded, Alphanumeric coded and 2D coded (Data-Matrix and TraXis). For different sizes in a 96-well configuration, adapter blocks are required for height adjustment. For the 0.50ml and 0.75ml tubes, we offer another type of Manual Capmat Sealer to guarantee easy capping.

The Micronic Capmats and Capclusters are recommended for the sealing of your valuable samples.

The Capmat exists out of a retaining film that holds 96 fixed TPE caps for sealing 96 tubes in a rack.

The Capcluster enables you to seal 96 tubes as easily as applying a Capmat, simply by removing the retaining film in order to leave each tube individually capped.

The functionality of the Capmat sealer is developed in conjunction with our users. Micronic has more than 25 years of experience in the life science market.

For more information:  
e-mail to [sales@micronic.com](mailto:sales@micronic.com)


**Micronic Europe B.V.**  
(Asia, Africa, Europe and Oceania)  
Phone: +31 (0)320 277 070  
**Micronic America LLC**  
(North, Central and South America)  
Phone: +1 484 480 3372

Please visit our website: [www.Micronic.com](http://www.Micronic.com)



#### INFORMATION

- Standard for all 1.4ml tubes
- Easy installation
- Integrity of your sample ensured
- Robust design
- One slide capping method




## MICRONIC | SAMPLE HANDLING | CAPPING & DECAPPING

### SPECIFICATIONS


<b>MATERIAL:</b>	<b>STAINLESS STEEL</b>
<b>DIMENSIONS:</b>	<b>350CM L x 206 H x 114 W</b>
<b>WEIGHT:</b>	<b>6 KG</b>

### TECHNICAL DRAWING

**MP55006**



**MP55005**



### INFORMATION

- Standard for all 1.4ml tubes
- Easy installation
- Integrity of your sample ensured
- Robust design
- One slide capping method

### ORDERING INFORMATION

Adapters need to be ordered separately for tubes different then 1.40ml

<b>Manual Capmat Sealer 0.50ml</b>		
MP55006	Manual Capmat Sealer	unit
MP55702	Manual Capmat Sealer Adapter for 0.50ml tubes	unit
<b>Manual Capmat Sealer 0.75ml</b>		
MP55006	Manual Capmat Sealer	unit
MP55700	Manual Capmat Sealer Adapter for 0.75ml tubes	unit
<b>Manual Capmat Sealer 1.10ml</b>		
MP55006	Manual Capmat Sealer	unit
MP55701	Manual Capmat Sealer Adapter for 1.10ml tubes	unit
<b>Manual Capmat Sealer 1.40ml</b>		
MP55005	Manual Capmat Sealer	unit



Micronic has a broader product range then the products listed under ordering information. For our complete product portfolio, please send your e-mail to [sales@micronic.com](mailto:sales@micronic.com) or contact us by phone.