

MICRONIC I SAMPLE HANDLING I EQUIPMENT & SOFTWARE

TUBE SELECTOR

IS IT DIFFICULT TO OBTAIN THE CORRECT SAMPLE FROM YOUR RACK? IF SO, THEN MICRONIC OFFERS YOU THE TUBE SELECTOR FOR 'HANDS-FREE' MANUAL PICKING AND PLACING OF SAMPLE TUBES FROM YOUR STORAGE RACK. COMPATIBLE WITH ALL STANDARD ANSI/SLAS FORMAT 96-RACK TUBES, THE VERSATILE TUBE SELECTOR REDUCES SAMPLE CONTAMINATION, AS CHERRY PICKING BY HAND IS NO LONGER NECESSARY.

Offering the convenience of singlehanded operation, the Tube Selector is very easy to use. A single thumb button enables the rapid and secure picking of a storage tube. Pressing the thumb button again releases the tube from the Tube Selector.

The Micronic tube range provides a storage solution for the samples stored in your lab. For secure tracking of your samples, the tubes are available with a laser etched 2D code (Data-Matrix or TraXis) on the bottom of the tube. This ensures traceability when using these tubes in your lab.

Micronic tubes in a 96-well format are compatible with the Micronic TPE or Screw Caps (TPE Caps also fit Screw Cap tubes). The tubes are available in bulk, refill carrier-96, Loborack-96 and Roborack-96.



Both racks have the 'Twist-Lock' feature to prevent turning during Screw Capping.

Micronic offers a range of low cost Starterpacks for organisations looking to implement traceable sample storage work practices. A Starterpacks contains everything needed to start using alphanumeric or 2D coded sample storage tubes.

Each Starterpack contains an easyto-use tube scanner, operating software, a case of racked tubes, (multi-piercable) tube caps, a versatile plate sealing device and manual tube decappers.

For more information: e-mail to sales@micronic.com

Micronic Europe B.V.

(Asia, Africa, Europe and Oceania)
Phone: +31 (0)320 277 070

Micronic America LLC

(North, Control and South America)

(North, Central and South America)


Phone: +1 484 480 3372

Please visit our website: www.Micronic.com



INFORMATION

- Manual thumb operated ejection
- Reduced contamination
- Single-handed operation
- Wipe clean surface
- Easy grip





MICRONIC I SAMPLE HANDLING I EQUIPMENT & SOFTWARE

SPECIFICATIONS

MATERIAL POLYPROPYLENE

HEIGHT 15CM DIAMETER 1.6CM



ORDERING INFORMATION

Tube Selector		
MP54005	Tube Selector	box/5
MP66000	Tube Selector	box/1



Micronic has a broader product range then the products listed under ordering information. For our complete product portfolio, please send your e-mail to sales@micronic.com or contact us by phone.